



Roger™

Hear better over noise and distance

A Sonova brand

PHONAK
life is on

Bridging the understanding gap

Whether you're in a restaurant, at work or relaxing at home, understanding speech is important to leading a fulfilled life.

Today's hearing aid Technology does an excellent job of improving speech understanding. However, even the most advanced hearing aids have limitations - especially in noisy places or at a distance. In such cases, state-of-the-art wireless microphones are required to boost performance and bridge the understanding gap. This is where the Roger microphone portfolio truly excels.

The Roger portfolio consists of microphones and receivers that function by picking up the voice of a speaker and wirelessly transmitting it to the listener while reducing background noise. This enables you to clearly hear speech in meetings or social gatherings, making hearing a pleasure.



Hear better in noisy situations with MultiBeam Technology


The theory is simple: the more spoken words you can understand, the more you will enjoy communication.

MultiBeam Technology is an innovative Roger microphone feature from Phonak. It enables you to experience exceptional speech understanding in group conversations such as business meetings or in noisy restaurants. By utilizing multiple microphones in six directions, speech from a 360° radius is calculated and compared. The direction with the best clarity is automatically selected. This results in a XX% improvement in speech intelligibility in a group conversation when background noise is present, compared to hearing aids alone*.

The MultiBeam Technology is now available in Roger Table Mic II and Roger Select.

* Speech perception increase by XX%* in a group conversation compared to using HI alone in noisy situations.

Linda Thibodeau study.



When you can focus on work rather than on hearing, life is on

Roger Table Mic II is a microphone dedicated for working adults who participate in various meetings. It automatically selects the person who's talking while reducing background noise and the effects of distance between you and the speaker.

Fully participate in large meetings

To cover a large meeting, several Roger Table Mic II can work together in a MultiTalker Network. With an outstanding range of more than 20 meters between you and the microphone, you will be able to fully participate and won't miss a thing.



Designed for ease of use

Roger Table Mic II is fully automatic - which makes it very easy to use. With a battery capacity of more than 16 hours*, you won't have to worry about the operating time.

The discreet design makes it blend in perfectly with the rest of your meeting equipment. For extra discreet usage you can choose to operate the table microphone with a remote control.

Choose the way you hear in large meetings

Some meetings have more background noise than others. The adjustable pick-up range allows you to adapt your listening comfort and choose the way you hear in various environments. Decide whether to have a broad pick-up range in areas with little background noise or a more focused pick-up range in areas where background noise is present.

* With 16 hour autonomy Roger Table Mic II last a full working day according to measurements.

When you feel at the heart of the conversation, life is on

Roger Select is a versatile microphone ideal for stationary situations where background noise is present such as restaurants.

Feel part of the whole conversation

When placed in the center of a table, it discreetly and automatically selects the person who is talking and seamlessly switches from one talker to another.

Be in control, decide whom to listen to

In situations where multiple conversations are taking place, with a single manual touch, select the person to hear and the conversation to focus on.



Enjoy precious moments

Your friends and family can simply clip it on their clothing or place it on a neck loop so you can clearly capture all those precious moments.

Experience your surroundings

Roger Select features wideband Bluetooth for phone calls and it can stream sound from TVs, music players, computers and more. Working with miniature Roger receivers, there is one for almost every hearing aids and cochlear implant sound processor.



Other Roger microphones

Besides the two new microphones: Roger Table Mic II and Roger Select, there are other advanced Roger microphones that offer a solution for your needs.



Roger Pen™

Roger Pen is a handy microphone for various listening situations. Thanks to its portable design, it can be conveniently used over distance and in noise. It has wideband Bluetooth for phone calls.



Roger EasyPen

Roger EasyPen is an effective yet simple microphone with only two buttons for clients who desire ultimate simplicity with full Roger performance.



Roger Clip-On Mic

Roger Clip-On Mic is a small microphone designed for one-to-one conversations. Featuring a directional microphone, the user can focus on a conversation with their partner.



Receivers

With three types of receivers available there is a Roger solution for you, regardless of the type of hearing aids, cochlear implants or Baha you currently wear.



Roger design-integrated receivers

These are designed to click onto Phonak hearing aids. They are also available for selected cochlear implants from Advanced Bionics, Cochlear and MED-EL.



Roger X

This miniature universal Roger receiver is compatible with many behind-the-ear hearing aids, cochlear implant sound processors and can be connected to many streamers.



Roger MyLink

This Roger receiver with a neck-loop works with any hearing aid or cochlear implant that features a telecoil.

Roger at work

Modern work life is highly demanding, especially on people's ability to communicate. Meetings often cause hearing difficulties for people with a hearing loss. Some of the main obstacles are due to side conversations, surrounding noise or the distance between you and someone at the other end of the table. Regardless of the situation, Roger microphones are always easy to use and enable you to fully communicate, participate and contribute at work.




Small meetings

Even small meetings can become overwhelming for people wearing hearing aids. Simply place one table microphone on the meeting table and you will hear your colleagues better.




Large meetings

If you are participating in a large meeting you might need more than one table microphone. You can use as many Roger Table Mic II as necessary. They transmit speech to your hearing aids from wherever the voice is coming.




Meetings with presenters

Presenters often stand at a distance from the meeting table. When the presenter wears a Roger microphone that works cohesively with the table microphone you can hear the speech clearly from both the presenter and the meeting participants.




One-on-one conversation in noise

When surrounded by distracting noise, the most effective approach is to hang a Roger microphone around your colleague's neck or clip it on the collar.



Multimedia

Technologies like video conferencing systems and Internet calls are frequently used at work. Roger microphones make hearing audio easy. Simply connect your Roger microphone, using the audio cable provided to bring these essential sounds directly to your ears.



Phone calls

Making phone calls is a real challenge for people affected by hearing loss. By connecting the Roger Pen or Roger Select to your phone/computer via Bluetooth or cable you can hear phone calls directly through your hearing aids or cochlear implant.

Roger in your private life

Many people with hearing loss would love to hear well and understand speech in busy social settings. Wouldn't it be nice if you no longer feel the urge to avoid these situations? Roger Technology puts the 'life' back into your social life – helping you to stay connected regardless of the background noise or listening environment.



Dinner with friends or family

Dinner tables and restaurants can often pose a real listening challenge due to their complex mix of chatter, background music, and clinking cutlery and glasses. A solution such as Roger Select cuts out these distracting noises and helps you to capture every word.



At parties

When music is playing and people are talking and laughing, it is easier to focus on a conversation with Roger because it delivers the conversation directly into your ears. Roger Pen is ideal for these listening situations. With its pen-like design it looks great and is easy to use.



At home

At home Roger reduces distracting background noise like kids playing, kitchen aids or the blare of a TV. This allows you to fully enjoy conversations with family members and friends.



In the car

Engine sounds, road noise and passengers facing away from you, make conversations challenging. By bringing the speaker's words directly into your ears, Roger makes it possible to overcome the communication challenges. Hearing in the car becomes enjoyable again.



Phone calls

Enjoy talking on the phone by connecting the Roger Pen or Roger Select to your phone or computer via Bluetooth or cable.



Watching TV

Roger microphones can be easily connected to your TV's audio output. This makes it possible for you to hear your favorite programs without having to increase the TV set's main volume.

Learn which Roger microphone suits your listening needs



Take the next step to find out which Roger system is best for you, talk to your hearing care professional or read more at www.phonakpro.com/Easyguide





Get connected

Roger enables you to effortlessly participate and contribute with your colleagues at work, friends and relatives. Roger makes certain that you hear everything in high quality sound allowing you to focus on the conversation without distraction.

Reimbursement

Roger systems needed for work are often, partly or fully, reimbursed. Sometimes reimbursement is also offered for solutions that cover your private life. Hearing care professionals offering the Roger portfolio can help you apply for reimbursement and support you in the application process.

Learn more about which Roger solution meets your needs. Talk to your hearing care professional or read more at www.phonak.com.



Join us online

There is a host of useful Web pages and tools online. Why not try them today?

HearingLikeMe

An online community for parents, children, teens and anyone whose life has been touched by hearing loss.

www.hearinglikeme.com

Facebook

The very latest Phonak updates in one handy, shareable place.

www.facebook.com/Phonak

Twitter

www.twitter.com/phonak

YouTube

The official Phonak YouTube channel is full of inspiring and educational content – from product movies and celebrity ambassadors to those heart-warming moments when an infant hears for the first time.

www.youtube.com/phonakofficial

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonak.com